


2022 order form

Use dimension chart to select appropriate size.



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

Conversion kit K116 \$745

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.	Important: User's weight must not exceed the maximum working load.			
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25 x 25	25 x 25½	32 x 29	39 x 33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175



* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color:
pink, blue or lime

		K110 mini	K130 small	K150 medium	K170 large
CHOOSE	K101	\$2,000 <input type="checkbox"/>	K103 \$2,260 <input type="checkbox"/>	K105 \$2,730 <input type="checkbox"/>	K107 \$3,260 <input type="checkbox"/>
	Pink	<input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>
	Blue	<input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>
	Lime	<input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>
Wheels (pair) (include brakes)		K113 \$1,200 <input type="checkbox"/>	K113 \$1,200 <input type="checkbox"/>	K115 \$1,355 <input type="checkbox"/>	K117 \$1,630 <input type="checkbox"/>
Mounting bracket		K133 \$68 <input type="checkbox"/>	K133 \$68 <input type="checkbox"/>	K173 \$79 <input type="checkbox"/>	K173 \$79 <input type="checkbox"/>
Arm prompts (pair) (require mounting bracket)		K512 \$510 <input type="checkbox"/>	K512 \$510 <input type="checkbox"/>	K512 \$510 <input type="checkbox"/>	K512 \$510 <input type="checkbox"/>
Hand loops (require mounting bracket)	Single	K578 \$84 <input type="checkbox"/>	K578 \$84 <input type="checkbox"/>	K578 \$84 <input type="checkbox"/>	K578 \$84 <input type="checkbox"/>
	Pair	K579 \$168 <input type="checkbox"/>	K579 \$168 <input type="checkbox"/>	K579 \$168 <input type="checkbox"/>	K579 \$168 <input type="checkbox"/>
Communication tray (requires mounting bracket)		K518 \$385 <input type="checkbox"/>	K518 \$385 <input type="checkbox"/>	K518 \$385 <input type="checkbox"/>	K518 \$385 <input type="checkbox"/>
Additional seat pad	Small	N/A	N/A	K135 \$147 <input type="checkbox"/>	K135 \$147 <input type="checkbox"/>
	Large	N/A	N/A	K155 \$158 <input type="checkbox"/>	K155 \$158 <input type="checkbox"/>
Sandals with wedges (pair)	Small	K114 \$110 <input type="checkbox"/>	K114 \$110 <input type="checkbox"/>	N/A	N/A
	Medium	N/A	N/A	K154 \$121 <input type="checkbox"/>	N/A
	Large	N/A	N/A	N/A	K174 \$126 <input type="checkbox"/>
Total		K110 \$ _____	K130 \$ _____	K150 \$ _____	K170 \$ _____

* Required