

# Choose your Anchor

Users with low tone receive support and develop strength for improved functional sitting. Users with high tone or movement disorders gain stability for increased functional use of the free hand and arm. Anchors promote inclusion and social development for the user.

**Anchors are short-term positioning aids and are not intended for primary standing support or transfers.**



**Hand Anchor:** This ergonomic hand Anchor attaches to any smooth surface to achieve stability and positioning. No more drilling through wheelchair trays!

K821 \$158  
Overall length: 5½"  
Length of handhold with bulb: 4¼"  
Diameter of handhold: 1¼"



**Grab Bar:** The horizontal grab bar attaches to any smooth, horizontal surface to assist in functional activities. It can be used to stabilize both hands at once.

K824 \$210  
Overall length: 20"  
Overall length with cups in: 13½"  
Diameter of handhold: 1¼"  
Length of bar: 10"



**Elbow Anchor:** The elbow Anchor provides stability to the forearm while leaving the hands free for functional activities. It can also reduce abnormal posturing of the shoulder girdle by maintaining the elbow position.

**Right** K815 \$184  
**Left** K816 \$184  
Overall length: 11"  
Overall width: 7¼"  
Width of arm trough: 4"  
Length of arm trough base: 9½"



**Wrist Anchor:** For a more secure handhold and to help maintain hand placement, the wrist Anchor includes a wrist strap and pad.

K822 \$179  
Overall length: 8"  
Length of handhold with bulb: 4"  
Diameter of handhold: 1¼"



**Horizontal Anchor:** The versatile horizontal Anchor is ideal for clients who can grasp more easily with their forearm pronated. Attach to any smooth table or tray surface for stability and positioning.

K823 \$168  
Overall length: 9½"  
Length of handhold with bulb: 4½"  
Diameter of handhold: 1¼"



**Arm Anchor:** The padded arm Anchor stabilizes the forearm. The handgrip extends and rotates to any angle to meet multiple positioning needs.

**Small right** K811 \$200  
**Small left** K812 \$200  
Inside of thumb to elbow: 7½"–12"  
Overall length with handhold: 9"–13"  
Width of arm trough: 4"  
**Large right** K813 \$230  
**Large left** K814 \$230  
Inside of thumb to elbow: 10"–15"  
Overall length with handhold: 11"–16"  
Width of arm trough: 4½"