


# 2023 order form

Use dimension chart to select appropriate size.



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

**Conversion kit** K116 \$785

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 <b>Key user dimension: height</b> Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.				
	<b>Important:</b> User's weight must not exceed the maximum working load.			
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25 x 25	25 x 25½	32 x 29	39 x 33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175



\* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

**Specify color:**  
pink, blue or lime

		K110 mini	K130 small	K150 medium	K170 large
<b>CHOOSE</b>		K101 \$2,100 <input type="checkbox"/>	K103 \$2,375 <input type="checkbox"/>	K105 \$2,870 <input type="checkbox"/>	K107 \$3,425 <input type="checkbox"/>
		Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>
		Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>
		Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>
<b>Wheels</b> (pair) (include brakes)		K113 \$1,260 <input type="checkbox"/>	K113 \$1,260 <input type="checkbox"/>	K115 \$1,425 <input type="checkbox"/>	K117 \$1,715 <input type="checkbox"/>
<b>Mounting bracket</b>		K133 \$70 <input type="checkbox"/>	K133 \$70 <input type="checkbox"/>	K173 \$85 <input type="checkbox"/>	K173 \$85 <input type="checkbox"/>
<b>Arm prompts</b> (pair) (require mounting bracket)		K512 \$535 <input type="checkbox"/>	K512 \$535 <input type="checkbox"/>	K512 \$535 <input type="checkbox"/>	K512 \$535 <input type="checkbox"/>
<b>Hand loops</b> (require mounting bracket)	Single	K578 \$90 <input type="checkbox"/>	K578 \$90 <input type="checkbox"/>	K578 \$90 <input type="checkbox"/>	K578 \$90 <input type="checkbox"/>
	Pair	K579 \$180 <input type="checkbox"/>	K579 \$180 <input type="checkbox"/>	K579 \$180 <input type="checkbox"/>	K579 \$180 <input type="checkbox"/>
<b>Communication tray</b> (requires mounting bracket)		K518 \$405 <input type="checkbox"/>	K518 \$405 <input type="checkbox"/>	K518 \$405 <input type="checkbox"/>	K518 \$405 <input type="checkbox"/>
<b>Additional seat pad</b>	Small	N/A	N/A	K135 \$155 <input type="checkbox"/>	K135 \$155 <input type="checkbox"/>
	Large	N/A	N/A	K155 \$165 <input type="checkbox"/>	K155 \$165 <input type="checkbox"/>
<b>Sandals with wedges</b> (pair)	Small	K114 \$115 <input type="checkbox"/>	K114 \$115 <input type="checkbox"/>	N/A	N/A
	Medium	N/A	N/A	K154 \$125 <input type="checkbox"/>	N/A
	Large	N/A	N/A	N/A	K174 \$135 <input type="checkbox"/>
<b>Total</b>		<b>K110</b> \$ _____	<b>K130</b> \$ _____	<b>K150</b> \$ _____	<b>K170</b> \$ _____

\* Required