


# 2023 order form

Use dimension chart to select appropriate size.



User dimensions (inches)	E420 small	E430 large
Height	30–50	46–72
 <b>Key user dimension: height</b> Select the appropriate stander by the user's overall height. Choose the model that allows for growth.		
	<b>Important:</b> User's weight must not exceed the maximum working load.	
Item dimensions (inches)	E420 small	E430 large
Board length and width	53 x 16½	73 x 19
Distance between laterals	7–13	9–16
Height when horizontal	28½	21–30
Height when vertical	54½	74
Base length and width	37½ x 27	46½ x 29½
Max. working load (lb)	100	250

\* **Supine Stander** includes padding, footboard, adjustable armrests, hip strap, chest strap, leg straps, one pair of laterals, one pair of head blocks and double-locking casters.

		E420 small	E430 large
		E40	\$3,130 <input type="checkbox"/>
		E41	\$4,900 <input type="checkbox"/>
<b>Additional laterals (pair)</b>	Small	E504	\$160 <input type="checkbox"/>
	Large	N/A	E439 \$235 <input type="checkbox"/>
<b>Tray</b>		E404	\$315 <input type="checkbox"/>
<b>CHOOSE</b>	<b>Adjustable abduction wedge</b>	E945	\$260 <input type="checkbox"/>
	<b>Round abductor</b>	E933	\$110 <input type="checkbox"/>
<b>Collar (requires round abductor)</b>		E934	\$90 <input type="checkbox"/>
<b>Sandals (pair)</b>	Small	E422	\$100 <input type="checkbox"/>
	Medium	E433	\$110 <input type="checkbox"/>
	Large	N/A	E432 \$110 <input type="checkbox"/>
<b>Wedges (pair, require sandals)</b>	Small	E424	\$70 <input type="checkbox"/>
	Medium	E435	\$75 <input type="checkbox"/>
	Large	N/A	E434 \$80 <input type="checkbox"/>
<b>Hand Anchor</b>		K821	\$165 <input type="checkbox"/>
<b>Total</b>		<b>E420</b>	<b>E430</b>
		\$ _____	\$ _____

\* Required